

with/without Pyridoxalphosphate

Diagnostic reagent for quantitative in vitro determination of ALAT (GPT) in serum or plasma on photometric systems

Order Information

Cat. No.	Kit size
01 00001 70 04 0125	R1 5 x 20 mL + R2 1 x 25 mL
01 00001 70 04 0500	R1 5x 80 mL + R2 1 x 100 mL
01 00001 70 10 0180	R1 4 x 36 mL + R2 4 x 9 mL
01 00001 70 02 0240	R1 4 x 36 mL + R2 4 x 9 mL
CDT-ALT	R1 3x 30 mL + R2 2x 11,3 mL

For determination with Pyridoxalphosphate activation additionally required:

01 00120 70 04 0018	6 x 3 mL
---------------------	----------



Summary [1,2]

Alanine Aminotransferase (ALAT/ALT), formerly called Glutamic Pyruvic Transaminase (GPT) and Aspartate Aminotransferase (ASAT/AST), formerly called Glutamic Oxalacetic Transaminase (GOT) are the most important representatives of a group of enzymes, the aminotransferases or transaminases, which catalyze the conversion of α -keto acids into amino acids by transfer of amino groups. As a liver specific enzyme, ALAT is only significantly elevated in hepatobiliary diseases. Increased ASAT levels, however, can occur in connection with damages of heart or skeletal muscle as well as of liver parenchyma. Parallel measurement of ALAT and ASAT is, therefore, applied to distinguish liver from heart or skeletal muscle damages. The ASAT/ALAT ratio is used for differential diagnosis in liver diseases. While ratios < 1 indicate mild liver damage, ratios > 1 are associated with severe, often chronic liver diseases.

Method

Optimized UV-test according to IFCC (International Federation of Clinical Chemistry and Laboratory Medicine)[modified].

Principle



Addition of pyridoxal-5-phosphate (P-5-P), recommended by IFCC stabilizes the activity of transaminases and avoids falsely low values in samples containing insufficient endogenous P-5-P, e.g. from patients with myocardial infarction, liver disease and intensive care patients [1, 3].

Reagents

Components and Concentrations

R1:	TRIS	pH 7.15	140 mmol/L
	L-Alanine		700 mmol/L
	LDH (lactate dehydrogenase)		≥ 2300 U/L
R2:	2-Oxoglutarate		85 mmol/L
	NADH		1 mmol/L

Pyridoxalphosphate

Good's buffer	pH 9.6	100 mmol/L
Pyridoxal-5-phosphate		13 mmol/L

Storage Instructions and Reagent Stability

The reagents are stable up to the end of the indicated month of expiry, if stored at 2 – 8°C, protected from light and contamination is avoided. Do not freeze the reagents!

Warnings and Precautions

- The reagents contain sodium azide (0.95 g/L) as preservative. Do not swallow! Avoid contact with skin and mucous membranes.
- Reagent 1 contains animal material. Handle the product as potentially infectious according to universal precautions and good clinical laboratory practices.
- In very rare cases, samples of patients with gammopathy might give falsified results [6].

- Sulfasalazine and sulfapyridine medication may lead to false results in patient samples. Blood collection must be done before drug administration.
- Please refer to the safety data sheets and take the necessary precautions for the use of laboratory reagents. For diagnostic purposes, the results should always be assessed with the patient's medical history, clinical examinations and other findings.
- For professional use only!

Waste Management

Please refer to local legal requirements.

Reagent Preparation

Substrate Start

The reagents are ready to use.

For the determination with Pyridoxalphosphate (P-5-P) mix 1 part of P-5-P with 100 parts of reagent 1, e.g. 100 μ L P-5-P + 10 mL R1

Stability after mixing:	6 days	at	2 - 8 °C
	24 hours	at	15 - 25 °C

Sample Start

without Pyridoxalphosphate

Mix 4 parts of R1 + 1 part of R2

(e.g. 20 mL R1 + 5 mL R2) = mono reagent

Stability:	4 weeks	at	2 - 8 °C
	5 days	at	15 - 25 °C

The mono reagent must be protected from light!

Materials required but not provided

DiaSystem Pyridoxalphosphate in case of determination with P-5-P activation (Cat.-no. 01 00120 70 04 0018)

NaCl solution 9 g/L

General laboratory equipment

Specimen

Serum, heparin plasma or EDTA plasma

Stability [4]:

3 days	at	20 - 25 °C
7 days	at	4 - 8 °C
7 days	at	-20 °C

Only freeze once! Discard contaminated specimens!

Assay Procedure

Application sheets for automated systems are available on request.

Wavelength	340 nm, Hg 365 nm, Hg 334 nm
Optical path	1 cm
Temperature	37°C
Measurement	Against air

Substrate Start

Sample or calibrator	100 μ L
Reagent 1	1000 μ L
Mix, incubate for 5 min., then add:	
Reagent 2	250 μ L
Mix, read absorbance after 1 min and start stopwatch.	
Read absorbance again 1, 2 and 3 min thereafter.	

Sample Start

Do not use sample start with Pyridoxalphosphate!

Sample or calibrator	100 μ L
Mono reagent	1000 μ L
Mix, read absorbance after 1 min. and start stopwatch.	
Read absorbance again 1, 2 and 3 min thereafter.	

Calculation

With factor

From absorbance readings calculate $\Delta A/\text{min}$ and multiply by the corresponding factor from table below:

$$\Delta A/\text{min} \times \text{factor} = \text{ALAT activity [U/L]}$$

	Substrate Start	Sample Start
340 nm	2143	1745
334 nm	2184	1780
365 nm	3971	3235

With calibrator

$$\text{ALAT [U/L]} = \frac{\Delta A/\text{min Sample}}{\Delta A/\text{min Calibrator}} \times \text{Conc. Calibrator [U/L]}$$

Conversion factor

$$\text{ALAT [U/L]} \times 0.0167 = \text{ALAT [\mu\text{kat/L}]}$$

Calibrators and Controls

For the calibration of automated photometric systems, DiaSystem UniCal CC calibrator is recommended. This method has been standardized against the original IFCC formulation. For internal quality control, Diasystem UniLab N and UniLab P controls should be assayed. Each laboratory should establish corrective action in case of deviations in control recovery.

	Cat. No.	Kit size
UniCal CC	06 00122 70 04 0018	6 x 3 mL
UniLab N	07 00123 70 05 0030	6 x 5 mL
UniLab P	07 00124 70 05 0030	6 x 5 mL

Performance Characteristics

Measuring range

On automated systems the test is suitable for the determination of ALAT activities up to 600 U/L. In case of a manual procedure, the test is suitable for ALAT activities which correspond to a maximum of $\Delta A/\text{min}$ of 0.16 at 340 and 334 nm or 0.08 at 365 nm. If such values are exceeded the samples should be diluted 1 + 9 with NaCl solution (9 g/L) and results multiplied by 10.

Specificity/Interferences

No interference was observed by ascorbic acid up to 30 mg/dL, bilirubin up to 40 mg/dL, hemoglobin up to 400 mg/dL and lipemia up to 2000 mg/dL triglycerides. For further information on interfering substances refer to Young DS [6].

Sensitivity/Limit of Detection

The lower limit of detection is 4 U/L.

Precision

Without Pyridoxal-5-phosphate

Intra-assay precision n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	22.2	1.38	6.22
Sample 2	44.8	1.17	2.62
Sample 3	101	1.02	1.00

Inter-assay precision n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	22.8	0.70	3.08
Sample 2	42.6	0.68	1.60
Sample 3	99.3	0.92	0.92

With Pyridoxal-5-phosphate

Intra-assay precision n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	33.8	1.25	3.71
Sample 2	72.0	2.04	2.83
Sample 3	128	2.77	2.16

Inter-assay precision n = 20	Mean U/L]	SD [U/L]	CV [%]
Sample 1	33.3	0.99	2.96
Sample 2	72.1	1.36	1.88
Sample 3	133	1.76	1.32

Method Comparison

With Pyridoxal-5-phosphate

A comparison of DiaSystem ALAT (GPT) FS with P-5-P (y) and the IFCC reference reagent (x) using 51 samples gave following results:
 $y = 1.000x - 0.200 \text{ U/L}$; $r = 0.999$

A comparison of DiaSystem ALAT (GPT) FS with P-5-P (y) and a commercially available test (x) using 51 samples gave following results: $y = 0.970x + 0.531 \text{ U/L}$; $r = 1.000$

Without Pyridoxal-5-phosphate

A comparison of DiaSystem ALT (IFCC) without Pyridoxalphosphate (y) with a commercially available test (x) using 51 samples gave following results:

$$y = 0.971x + 0.047 \text{ U/L}; r = 1.000$$

Reference Range

With pyridoxal-5-phosphate activation

Women [7]	< 34 U/L	< 0.57 $\mu\text{kat/L}$
Men [7]	< 45 U/L	< 0.75 $\mu\text{kat/L}$
Children [1]	1 – 30 day(s)	< 25 U/L < 0.42 $\mu\text{kat/L}$
	2 – 12 months	< 35 U/L < 0.58 $\mu\text{kat/L}$
	1 – 3 year(s)	< 30 U/L < 0.50 $\mu\text{kat/L}$
	4 – 6 years	< 25 U/L < 0.42 $\mu\text{kat/L}$
	7 – 9 years	< 25 U/L < 0.42 $\mu\text{kat/L}$
	10 – 18 years	< 30 U/L < 0.50 $\mu\text{kat/L}$

Without pyridoxal-5-phosphate activation

Women [8,9]	< 31 U/L	< 0.52 $\mu\text{kat/L}$
Men [8,9]	< 41 U/L	< 0.68 $\mu\text{kat/L}$

Each laboratory should check if the reference ranges are transferable to its own patient population and determine own reference ranges if necessary.

Literature

1. Thomas L. Alanine aminotransferase (ALT), Aspartate aminotransferase (AST). In: Thomas L, editor. Clinical Laboratory Diagnostics. 1st ed. Frankfurt: TH-Books Verlagsgesellschaft; 1998. p. 55-65.
2. Moss DW, Henderson AR. Clinical enzymology. In: Burtis CA, Ashwood ER, editors. Tietz Textbook of Clinical Chemistry. 3rd ed. Philadelphia: W.B Saunders Company; 1999. p. 617-721.
3. Bergmeyer HU, Horder M, Rej R. Approved Recommendation (1985) on IFCC Methods for the Measurement of Catalytic Concentration of Enzymes. L.Clin. Chem. Clin. Biochem 1986; 24: 481-495.
4. Bakker AJ, Mücke M. Gammopathy interference in clinical chemistry assays: mechanisms, detection and prevention. ClinChemLabMed 2007;45(9):1240-1243.
5. Guder WG, Zawta B et al. The Quality of Diagnostic Samples. 1st ed. Darmstadt: GIT Verlag; 2001; p. 14-5.
6. Young DS. Effects of Drugs on Clinical Laboratory Tests. 5th ed. Volume 1 and 2. Washington, DC: The American Association for Clinical Chemistry Press 2000.
7. Schumann G, Bonora R, Ceriotti F, Féraud G et al. IFCC primary reference procedure for the measurement of catalytic activity concentrations of enzymes at 37 °C. Part 4: Reference procedure for the measurement of catalytic concentration of alanine aminotransferase. Clin Chem Lab Med 2002;40: 718-24.
8. Lorentz K, Röhle G, Siekmann L. Einführung der neuen Standardmethoden 1994 zur Bestimmung der katalytischen Enzymkonzentrationen bei 37 °C. DG Klinische Chemie Mitteilungen 26; 1995; Heft 4.
9. Zawta B, Klein G, Bablok W. Temperature Conversion in Clinical Enzymology? Klin. Lab. 1994; 40: 33-42.

Manufacturer

DiaSystem Scandinavia AB
 Datorgatan 3, Sweden – 561 33 Jönköping
 Phone +46 36 126220 • Fax +46 36 187730
 info@diasystem.se • www.diasystem.se

