CANINE & FELINE T4

For veterinary use only

INTENDED USE

The Vcheck T4 is an *in vitro* diagnostic test kit for the quantitative measurement of T4 (thyroxine) concentration in canine serum or feline serum. T4 is one of the hormones secreted by the thyroid gland that is responsible for the regulation of metabolism. The Vcheck T4 can be used to assist in the monitoring of structural or functional abnormalities of the thyroid gland. The Vcheck T4 is designed to be used only by veterinarians.

PRINCIPLE

The Vcheck T4 Test Kit is based on the competitive immunoassay method for the quantitative measurement of T4 concentration. When the sample is delivered to an assay diluent tube, the sample and a fixed amount of colloidal gold-labeled anti-T4 monoclonal antibody are mixed in the assay diluent buffer using the disposable tablet pipette. T4 in the sample reacts with the gold-labeled antibody. Any unbound gold-labeled antibody in the mixture binds to T4-BSA coated on the membrane. The density of the test line is inversely proportional to the T4 concentration of the sample. The Vcheck Analyzer reads the density of the test line and calculates the T4 concentration from the calibration curve data. The control line is a reference line that indicates the test has been performed correctly.

MATERIALS PROVIDED

Reagent	10 Tests/Kit
① Vcheck T4 Test device	10
② Assay diluent tube	10
③ Disposable pipette tip	20
④ Disposable tablet pipette (yellow bulb)	10
⑤ Instructions for use	1

MATERIALS REQUIRED, BUT NOT PROVIDED

- BIONOTE Vcheck Analyzer
- 50 µl pipette
- 100 µl pipette

STORAGE AND STABILITY

- Store the test kit at 2~8°C. **DO NOT FREEZE.**
- Do not store the test kit in the direct sunlight.
- The test kit is stable until the expiry date that is marked on the package label.

PRECAUTIONS

- The test kit is for canine and feline use only. Do not use for other animals.
- The test device is sensitive to humidity and heat. Perform the test immediately after removing the test device from the aluminum foil pouch.
- Do not reuse the test components.
- Do not touch the membrane in the result window of the test device.
- Do not use the test kit beyond the stated expiry date marked on the label.
- Do not use the test kit if the pouch is damaged or the seal is broken.
- Do not mix components from different lot numbers; the components in this kit have been quality control tested as a standard batch unit.
- All samples should be handled as being potentially infectious. Wear protective gloves while handling samples. Wash hands thoroughly afterwards.
- Decontaminate and dispose of all samples, used kits, and potentially contaminated materials safely in accordance with national and local regulations.

- It is recommended to use a plain tube for sample collection. Test results cannot be guaranteed for anything other than the recommended tubes.

COLLECTION AND PREPARATION OF SAMPLE

- Canine serum or feline serum should be used with this test. A method for preparing the serum is as follows.
- [Serum]** Collect the whole blood using a plain tube **ONLY**. Allow to settle for 30 minutes for blood coagulation and then centrifuge to obtain a serum supernatant.
- If serum samples are not tested immediately, they should be refrigerated at 2~8°C and used within 7 days. For longer storage, serum can be frozen (-20°C or colder). Frozen samples should be brought to room temperature (15~30°C) prior to use.

TEST PROCEDURE

- Allow all kit components and sample to reach room temperature (15~30°C) prior to testing.
- Test result can be measured by two available testing procedures, "Incubate and Read" or "Read Immediately". Users can select one of the test procedures depending on the situation.

Incubate and Read


- Turn on V200 Analyzer and select "Standard Test".

- Remove the test device from the aluminum foil pouch. Once the "Insert Device" is displayed in the screen, insert the test device.

- Select "Canine" or "Feline" depending on the type of sample used.

- "Incubate and Read" and "Read Immediately" will appear on the display of the V200 Analyzer.


- Use a 50 µl pipette to draw 50 µl of the sample and add to the assay diluent tube (125 µl).

- Use a disposable tablet pipette (yellow bulb) to mix and dissolve the tablet in the assay diluent tube. Mix over 8 times until the tablet is completely dissolved. While mixing the tablet, squeeze the yellow bulb gently to minimize foam.

- Select "Incubate & Read" and wait 10 minutes for incubation. 10 seconds before the end of the incubation phase, an alarm will sound. \*Sample mixture should be incubated at room temperature (15~30°C).

- Draw 100 µl of diluted sample using a 100 µl pipette.

- Add 100 µl of the mixture into the sample hole of the test device and press START to initiate testing.


- The V200 analyzer will display the test result on the screen after another 10 minutes.

- Remove the test device.

Read Immediately

- Use a 50 µl pipette to draw 50 µl of the sample and add to the assay diluent tube (125 µl).

- Use a disposable tablet pipette (yellow bulb) to mix and dissolve the tablet in the assay diluent tube. Mix over 8 times until the tablet is completely dissolved. While mixing the tablet, squeeze the yellow bulb gently to minimize foam.

- Incubate the sample mixture at room temperature (15~30°C) for 10 minutes.


- During incubation of the sample, turn on V200 Analyzer and select "Standard Test".


- Remove the test device from the aluminum foil pouch. Once the "Insert Device" is displayed in the screen, insert the test device.

- Select "Canine" or "Feline" depending on the type of sample used.

- "Incubate and Read" and "Read Immediately" will appear on the display of the V200 Analyzer.

- Select "Read Immediately".

- After incubating the sample mixture for 10 minutes, draw 100 µl of diluted sample using a 100 µl pipette.

- Add 100 µl of the mixture into the sample hole of the test device and press START to initiate testing.

- The V200 analyzer will display the test result on the screen after another 10 minutes.

- Remove the test device.

INTERPRETATION OF THE RESULT

- Read the concentration value of canine T4 or feline T4 appearing on the display of the BIONOTE Vcheck Analyzer. (0.5~8.0 µg/dL)
- If "< 0.5 µg/dL" appears on the display, it means the concentration of canine T4 or feline T4 in the specimen is less than 0.5 µg/dL.
- If "> 8.0 µg/dL" appears on the display, it means the concentration of canine T4 or feline T4 in the specimen is greater than 8.0 µg/dL.
- 1 µg/dL is equal to 12.87 nmol/L.

\* It is recommended that each laboratory establishes its own reference range.

[Canine]

< 1.0 µg/dL	1.0~2.0 µg/dL	1.0~4.0 µg/dL	> 4.0 µg/dL
Low	Low normal	Normal	High

[Feline]

< 0.8 µg/dL	0.8~4.7 µg/dL	2.3~4.7 µg/dL	> 4.7 µg/dL
Subnormal	Normal	Normal (Gray zone)	High (hyperthyroidism)

SCREEN MESSAGES AND TROUBLE SHOOTING

[V200]

Error message	Error description
Contaminated Device	The test device is damaged or inserted improperly. <b>Solution:</b> Discard the test device and retest with a new test device and a new specimen.
Insufficient Sample	An insufficient amount of blood has been applied. <b>Solution:</b> Retest with a new test device with enough specimen, ensuring that blood is placed in to the narrow channel in the top edge of the test device.
Expired Device	The test devices are expired. <b>Solution:</b> Retest with a new test device that is not expired.
Temperature Error	The environmental temperature is above or below the operating range of the analyzer. <b>Solution:</b> Move to an area in the acceptable temperature range for the analyzer and perform the test. Do not heat or cool the analyzer artificially.
Printer Connection Fail	The communication between analyzer and barcode or printer has failed. <b>Solution:</b> Reconnect the analyzer and external device. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Barcode Error	The measured total hemoglobin is out of the range of 7 to 23 g/dL. <b>Solution:</b> This error occurs when a specimen has a total hemoglobin in the abnormal range. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Extremely High Total Hemoglobin	The test is invalid. <b>Solution:</b> Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Result: Invalid	The calibration is overdue. <b>Solution:</b> Retest with a new test device and a new patient specimen. If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
Calibration Overdue	A test device that is not supported by the analyzer has been loaded. <b>Solution:</b> Check whether the test device is manufactured by BioNote, Inc.
Not Supported Device	Internal error has occurred. <b>Solution:</b> If the error continues after turning ON/OFF the analyzer, please contact BioNote, Inc.
EEE	

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