VC.3023-VC.3024-HDS

camera





Contents

1.0 The Application of the VC.3023 and VC.3024-HDS Series Camera	3
2.0 Available ports on the back of the camera body (same for both cameras)	3
3.0 VC.3023 and VC.3024-HDS Series Camera Function Description	4
3.1 Video output	4
3.2 Image capture and video saving in SD card	4
3.3 Image adjustment function	4
3.4 Image operation function	4
3.5 Other functions	
4.0 VC.3023 and VC.3024-HDS series camera application configurations	5
4.1 Camera working standalone with built-in on-screen ImageFocus software	5
5.0 Brief introduction of VC.3023 and VC.3024-HDS series camera's UI and its functions_	6
5.1 On-screen ImageFocus user interface (UI)	6
5.2 The camera control panel on the Left Side of the Screen	7
5.3 The measurement toolbar on top of the screen	8
5.4 Icons and functions of the camera control toolbar at the bottom of the screen	9
6.0 Calibration process of HDMI cameras	14
6.1 Step 1	14
6.2 Step 2	15
6.3 Step 3	15
6.4 Step 4	
6.5 Step 5	16
7.0 Sample photos cantured with VC 3023 and VC 3024-HDS series camera	16

1.0 The application of the VC.3023 and VC.3024-HDS Series Camera

The VC.3023 and VC.3024-HDS series camera is intended to be used for the acquisition of digital images from the stereo microscope and biological microscope. The basic characteristics are listed below:

- · Sony Starvis back illuminated CMOS sensor
- · Full HD HDMI video outputs
- · SD card for captured images and video storage
- Embedded on-screen ImageFocus for the control of the camera
- With strong 'image adjustment' and other related processing functions



Figure 1.0 VC.3023 and VC.3024-HDS series camera

2.0 Available ports on the back of the camera body

(the same for both cameras)

Interface	Function Description
1. HDMI	Complies with HDMI 1.4 standard. 1080 P format video output for standard Full HD monitor
2. USB Mouse	Connect USB mouse for easy operation with embedded On-screen ImageFocus software
3. SD	Complies with SDIO 3.0 standard and SD card can be inserted for video and images storage
4. DC 12 V	Power adapter connection (12 V / 1 A)
5. LED	LED status indicator



Figure 2.0b VC.3024-HDS

3.0 VC.3023 and VC.3024-HDS series camera function description

3.1 Video output

Video output interface	Function description	
HDMI interface	Complies with HDMI 1.4 standard; 60fps@1080P	

3.2 Image capture and video saving in SD card

Function name	Function description
Video saving	Video format: 2 MP (1920 \times 1080) H264 encoded MP4 file: Video saving frame rate: $50\sim60$ fps (related with SD card performance);
Image capture	2 MP (1920 x 1080) JPEG image in SD card
Measurement saving	Measurement information saved in different layer with image content; Measurement information is saved together with image content in burn- in mode

3.3 Image adjustment function

Function name	Function description
Exposure / Gain	Automatic / Manual exposure
White balance	Manual / Automatic / ROI mode
Sharpening	Supported
3D denoise	Supported
Saturation adjustment	Supported
Contrast adjustment	Supported
Brightness adjustment	Supported
Gamma adjustment	Supported
50 HZ / 60 HZ anti-flicker function	Supported

3.4 Image operation function

Function name	Function description
Zoom In / Zoom Out	Up to 10X
Mirror / Flip	Supported
Freeze	Supported
Cross line	Supported
Embedded files browser	Supported
Video playback	Supported
Measurement function	Supported

3.5 Other functions

Function name	Function description
Restore factory settings	Supported
Multiple language support	English / Simplified Chinese / Traditional Chinese / Korean / Thai / French /
25-25 100 30 02-20	German / Japanese / Italian / Russian

4.0 VC.3023 and VC.3024-HDS series camera application configurations

4.1 Camera working standalone with built-in on-screen ImageFocus software

For this application, apart from the microscope, you only need an HDMI monitor, the camera comes with a USB mouse and the camera embedded on-screen ImageFocus software. The steps to start the camera are listed as below:

Connect the camera to a HDMI monitor using the HDMI cable (1)












Figure 4.1a VC.3023 and VC.3024-HDS series




Figure 4.1b Menus

5.0 Brief introduction of the camera's UI and its functions

5.1 On-screen ImageFocus user interface (UI)

The VC.3023 and VC.3024-HDS series camera's UI shown in Figure 5.1 includes a camera control panel on the left of the screen, a measurement toolbar on the top of the screen and a camera control toolbar on the bottom of the screen




Figure 5.1 The VC.3023 and VC.3024-HDS series camera control UI

Not	Notes		
1	To show the camera control panel, move your mouse to the left of the screen. See 4.1b and 5.1 for details		
2	Move the mouse cursor to the top of the screen, a measurement toolbar will pop up for calibration and measurement operations. When user left-clicks the Float/Fixed button on the measurement toolbar, this toolbar will be fixed. In this case the camera control panel will not pop up automatically, even if users move the mouse cursor to the left side of the screen. Only when user left-clicks the button on the measurement toolbar to exit from measuring procedure will he or she be able to do other operations on the camera control panel, or the camera control toolbar. During the measuring process, when a specific measuring object is selected, an object location & attributes control bar will appear for changing location and properties of the selected object. See Sec. 7.3 for details		
3	When users move mouse cursor to the bottom of the screen, the camera control toolbar will pop up automatically \bigcirc		

5.2 The camera control panel on the Left Side of the Screen

The camera control panel controls the camera to achieve the best video or image quality according to the specific applications; It will pop up automatically when the mouse cursor is moved to the left side of the screen. Left-clicking button to achieve Display / Auto hide switch of the camera control panel

amera control panel	Function	Function Description
3	Snap	Capture image and save it to the SD card
Camera Control Panel	Record	Record video and save it to the SD card
Snap Record	Auto exposure	When Auto exposure is checked, the system will automatically adjust exposure time and gain according to the value of exposure compensation
	Exposure compensation	Available when Auto exposure is checked. Slide to left or right to adjust Exposure compensation according to the current video brightness to achieve proper brightness value
Gain: 0	Exposure time	Available when Auto exposure is not checked. Slide to left or right to reduce or increase exposure time, adjusting brightness of the video
Red: 101	Gain	Adjust Gain to reduce or increase brightness of video The noise will be reduced or increased accordingly
Blue: 75	Red	Slide to left or right to decrease or increase the proportion of red in RGB on video
● Auto	Green	Slide to left or right to decrease or increase the proportion of green in RGB on video
Sharpness: 0	Blue	Slide to left or right to decrease or increase the proportion of blue in RGB on the video
Denoise: 0 Saturation: 50	Auto white balance	White balance adjustment according to the video continuously
Saturation: 50 Gamma: 6	Manual white balance	Adjust the red or blue item to set the video white balance
Contrast: 60	ROI white balance	White balance could be adjusted when the ROI regio is changed according to content inside the ROI regio
Brightness: 50	Sharpness	Adjust sharpness level of the video
O DC O AC(50Hz) AC(60Hz)	Denoise	Slide left or right to denoise the video
	Saturation	Adjust saturation level of the video
Default	Gamma	Adjust gamma level of the video. Slide to the right side to increase gamma and to the left to decrease gamma
	Contrast	Adjust contrast level of the video. Slide to the right side to increase contrast and to the left to decrease contrast
	DC	For DC illumination, there will be no fluctuation in light source, so no need for compensating light flickering
	AC (50HZ)	Check AC (50HZ) to eliminate flickering caused by 50Hz light source
	AC (60HZ)	Check AC (60HZ) to eliminate flickering caused by 60Hz light source
	Default	Restore all the settings in the camera control panel to default values

5.3 The measurement toolbar on top of the screen

The **measurement toolbar** will pop up when moving the cursor to any place near the upper edge of the screen. Here is the introduction of the various functions on the **measurement toolbar**:



Figure 5.3 The measurement toolbar on the upper side of the screen

Icon	Function
1	Float / fix switch of the Measurement Toolbar
<u></u>	
✓ Visible	Show / hide Measurement
VISIBIC	Objects
Nanometer (nm) 🕶	Select the desired
Nanometer (IIII)	measurement unit
	Select magnification
4X ~	for measurement after
	calibration
	Object select
	Angle
4	Aligie
	
/\	4-Points angle
/ \	
	Point
•	**************************************
	A 1 % 15
	Arbitrary line
>	3-Points line
/	
	Horizontal line
	Tiorizontal line
i e	Vertical line
l,	
\/	3-Points vertical line
X	
	D
//	Parallel
~	
	Rectangle
_	Ellipse
	Linpse
\bigcirc	5-Points ellipse
\mathbf{O}	

Icon	Function
Θ	Circle
0	3-Points circle
<u></u>	Annulus
P	Two circles and its center distance
g g	3-Points two circles and its center distance
\circ	Arc
Î	Text
\Diamond	Polygon
5	Curve
um	Scale bar
7	Arrow
8	Execute calibration to determine the corresponding relation between magnification and resolution, which will establish the corresponding relationship between measurement unit and the sensor pixel size. Calibration needs to be done with the help of a micrometer
export -	Export the measurement information to CSV file(*.csv)
B	Measurement setup

Icon	Function
B	Delete all the measurement objects
X	Exit from measurement mode

Icon	Function
< >AV ♣ 1111	When the measurement ends, left-click on a single measuring object and the object location & properties control bar will show up. The user can move the object by dragging the object with the mouse. But more accurate movement can be done with the control bar. The icons on the control bar mean move left, move right, move up, move down, color adjustment and delete

Note:

- When the user left-clicks the Display/Hide button 🖋 on the measurement toolbar, this toolbar will be fixed. In this case the camera control panel will not pop up automatically even if moving the mouse cursor to the left edge of the screen. Only when user left-clicks the X button on the measurement toolbar to exit from the measurement mode will he or she be able to do other operations with the camera control panel or the camera control toolbar
- When a specific measurement object is selected during the measurement process, the object location & attributes control bar < > A will appear for changing the object location and properties of the selected objects

5.4 Icons and functions of the camera control toolbar at the bottom of the screen



Icon	Function
+	Zoom in
\ominus	Zoom out
<u> </u>	Horizontal flip
	Vertical flip
(c-G	Color / gray

Icon	Function
•	Video freeze
#	Display cross line
	Browse images and videos on the SD card
X	Settings
(i)	Check the version of on- screen ImageFocus



The setting is relatively more complicated than the other functions. Here is more information about it:

5.4.1 Setting>Measurement




Figure 5.4.1 The measurement setup

Global	Precision	Used to set the number of digits after the decimal point of the measurement result	
Calibration	Line Width	Used for defining the width of the lines for calibration	
	Color	Used for defining the color of the lines for calibration	
	EndPoint type	Used for defining shape of the endpoints of lines for calibration: null means no endpoints, rectangle means rectangle type of endpoints. It makes alignment easier	
Point, Angle, Lin	e, Horizontal line, Vert	ical line, Rectangle, Circle, Ellipse, Annulus, Two circles, Polygon, Curve	
		along with the measurement command mentioned above will unfold the ibute settings to set the individual property of the measurement objects	

5.4.2 Setting>Magnification




Figure 5.4.2 Comprehensive magnification calibration settings page



Name	The name given to the current calibrated magnification. This typically is the objective magnification, like 4x or 10x etc. Also, other user-defined information could be added into the name. For example, microscope model, operator name, etc.
Resolution	Pixels per meter. Image device like microscopes have high resolution value
Clear All	Click the Clear All button to clear the calibrated magnifications
Delete	Click Delete to delete the selected magnification

5.4.3 Settings>Image Format




Figure 5.4.3 Comprehensive image format settings page

Measurement object save method	Burn in mode: The measurement objects are merged into the current image. User cannot edit the measurement objects anymore. This mode is irreversable
	Layered mode: The measurement objects are saved in a different layer with current image data in the target file. User can edit the measurement objects in the target file with photo editing software (not included) on the PC. This mode is reversable

5.4.4 Settings>Video




Figure 5.4.4a Comprehensive setting of video settings page-playback




Figure 5.4.4b Comprehensive setting of video settings page-video encode

Fast forward/Reverse interval	The time interval of the playback of video files	
Video encode	H264: The encoding format of the video files is H264 format	
	H265: The encoding format of the video files is H265 format	

5.4.5 Setting>Storage




Figure 5.4.5 Comprehensive setting of SD card setting page

Storage device	SD card: SD card is only supported as storage device	
File system format of the storage	List the file system format of the current storage device	
device	FAT32:	The file system of SD card is FAT32. The maximum video file size of single file is 4 GB
	exFAT:	The file system of SD card is exFAT. The maximum video file size of single file is 4 GB
	NTFS:	The file system of SD card is NTFS. The maximum video file size of single file is 4 GB.
		Use PC to format the SD cards and switch between FAT32, exFAT and NTFS
	Unknown St	tatus: SD card not detected or the file system is not identified

5.4.6 Setting>Files




Figure 5.4.6 Comprehensive setting of files settings page

Image File Name	Auto: The image files will be saved automatically with the specified prefix Manual: Users has to specify the file name before image saving
Video File Name	Auto: The video file will be saved automatically with the specified prefix Manual: Users has to specify the video file name before video recording
<i>Note</i> : The maximum vio	leo file size is 4 GB. Multiple video files may be generated automatically during long time video

5.4.7 Setting>Language



Figure 5.4.7 Comprehensive setting of language selection setting page

English	Set language of the whole software into English;
Simplified Chinese	Set language of the whole software into Simplified Chinese;
Traditional Chinese	Set language of the whole software into Traditional Chinese;
Korean:	Set language of the whole software into Korean;
Thailand	Set language of the whole software into Thailand;
French	Set language of the whole software into French
German	Set language of the whole software into German
Japanese	Set language of the whole software into Japanese
Italian	Set language of the whole software into Italian
Russian	Set language of the whole software into Russian

5.4.8 Setting>Miscellaneous





Figure 5.4.8 Comprehensive miscellaneous settings page

Auto exposure	The maximum exposure time during auto exposure process could be specified. Setting
	this item to a lower value could guarantee a faster frame rate during auto exposure
ROI Color	Choosing the ROI rectangle line color
Camera parameters import	Import the camera parameters from the SD card to use the previously exported camera parameters
Camera parameters export	Export the camera parameters to the SD card to use the previously exported camera parameters
Reset to factory defaults	Restore camera parameters to its factory status

6.0 Calibration process of HDMI cameras

6.1 Step 1

You need to put a calibration slide like below on your microscope stage:



6.2 Step 2

Click the calibration button which is circled out in red to start the calibration process



6.3 Step 3

After you click the **calibration button** a calibration window and a soft-keyboard will pop up. Use the soft-keyboard to enter magnification of the microscope and actual length of the calibration slide:

For example: when using a Euromex microscope at 4x: you enter in 4x.

Actual length is the actual length of the calibration slide you are using



You can also name the magnification to be "euromex_4x"



Click OK to complete the calibration process

6.4 Step 4

After the calibration process has been completed you will be able to find the magnification (4x) in menu bar:

\$\text{\$\frac{\pi}{2}\$ is visible Hillineter (nn) \$\sqrt{\pi}\$ \tag{\$\frac{\pi}{2}\$}\$ \$\frac{\pi}{2}\$ \$\frac

Before:

After:

6.5 Step 5

Change the unit of length/area in the dropdown list and start the measuring procedure



7.0 Sample photos captured with these cameras




Figure 7.0a Alfalfa stem



Figure 7.0b Top bud

Euromex Microscopen bv • Papenkamp 20 • 6836 BD Arnhem • The Netherlands T +31 (0) 26 323 22 11 • info@euromex.com • www.euromex.com





