# Evaluation of the Anigen Rapid Giardia Ag Test Kit

I submit this report to BioNote, Inc.

Yun young-min

College of Veterinary Medicine Jeju National University

### **■** Cross reaction test against intestinal parasites and bacteria

### **1. Kit**

[Anigen Rapid Giardia Ag Test Kit] Lot No. XXXXX

### 2. Intestinal parasites

- 1) Giadia lamblia, Coccidium
  - (It was kept in college of veterinary medicine, Jeju national university.)
- Toxocara sp., Anchylostoma sp.
   (It was extracted from infected dogs at the animal hospital of veterinary medicine, Jeju national university.)
- 3) E. coli, Salmonella sp.
- (It was isolated and identified from the bacteriological lab of veterinary medicine, Jeju national university.)


### 3. Method

- 1) Counted the number of cysts and bacteria. (When they were cysts, we diluted them with 10% suspension with PBS.)
- 2) Diluted them with assay diluents as same quantity.
- 3) Dropped them to the dropping site in the kit.
- 4) Observed the result within 10 minutes.

### 4. Result

All of the parasites (Coccidium) and bacteria were negative except the Giardia lamblia.

Name	Picture	Result
Giardia lamblia	Giardia Ag Giardia Giardia Giardia	Positive
Toxocara canis	Giardia Ag  toxocara  canis  09.2.12	Negative
Coccidium	Giardia Ag Cocc:dium 09. 2. 12	Negative
Anchylostoma canium	Giardia A9 AnchyRostona	Negative



## 5. Conclusion

[Anigen Rapid Giardia Ag Test Kit] reacted to Giardia lamblia only.

## ■ Detection limit according to the number of cysts

### **1. Kit**


[Anigen Rapid Giardia Ag Test Kit] Lot No. XXXXX

### 2. Giardia parasites

Giardia lamblia, kept in veterinary college, Jeju University, was inoculated to two dogs which were 60days old. It induced them to be infected with Giardia lamblia. We practiced Giardia lamblia test once a week and also did sensitivity of the kit according to the number of Giardia cysts.

### 3. Method

- 1) Floated Giardia cysts with assay diluents.
- 2) Counted them by the hemochamber.
- 3) Performed serial dilution. (62.5, 125, 250, 500, 1000, 2000, 4000, 8000cyst/100ul.)
- 4) Dropped 100ul to each kit.
- 5) Observed band density after 10 minutes later.
- 6) Measured correlation.



A: *Giardia* cyst in dilution of the kit B: *Giardia* cyst stained by iodine C,D: Giardia cyst in hemochamber (X400)

### 4. Result

Concentration	Band density	Picture
62.5/100ul	-	Giardia A9 62.5/10ml Gr. cycl
125/100ul	1+, 10%	Giardia A9  125/100ml  Gr. cyst
250/100ul	2+, 20%	Giardia Ag 250/coul Cr. cyst
500/100ul	2+, 30%	Giardia A9  500/100ul  Cr. cyst
1000/100ul	2+, 40%	Giardia Ag C T  /000/1004L  G. cyst
2000/100ul	2+, 50%	Giardia A9 2000/1001  Gr. cyst
4000/100ul	3+, 65%	Giardia Ag C T 4000/1000l  Gr. cysk
8000/100ul	3+, 80%	Giardia Ag  Seco/local  Cr. cyst

# 5. Conclusion

The detection limit of [Anigen Rapid Giardia Ag Test Kit] is 125 cysts/100 ul.

### 12. Comparison test with other commercial kit

### 1) Material

### (1) Using reagent

A. Testing reagent : Anigen Rapid Giardia Ag Test kit (Korea)

Lot No.	Manufactured date	Expiration date
T 8002	2008. 4. 23.	

### B. Compared reagent

Name	Manufacturing firm	Lot No.	
SNAP Giardia	IDEXX	09511-ED416	

### C. Microscopic Examination

The fresh feces were collected from the rectum. We coated them on the slide glass slightly and observed them through the microscope at x100 or x400.

#### (2) Samples

Total 176 examined that suspected as diarrhea feces from five animal hospitals.

### 2) Method

Anigen Rapid Giardia Ag Test kit was used for diarrhea feces from five animal hospitals. Microscopic examination or IDEXX kit test was regarded as the gold standard at the same time. And we calculated specificity and sensitivity as following methods.

3) Examiner: Five animal hospitals (Reference, 12. Refer to confirmation of field test)

#### 4) Result

Anigen Rapid Giardia Ag Test kit showed 100% of sensitivity and specificity. However, sensitivity and specificity of IDEXX are 95~96% and 99~100%. Therefore actual sensitivity and specificity of Anigen Rapid Giardia Ag Test kit are expected as over 95% and 99% (Lieb MS and Zajac AM. Giardiasis in dogs and cats. Veterinary Medicine. 1999:793-802, Hill SL *et al.* Prevalence of enteric zoonotic organ in cats. JAVMA. 2000;216(5):687-692).

Animal hospital name	Method	Total number examined	Difference test		Anigen Rapid Giardia Ag Test kit	
Hae-ma-ru animal hospital	IDEXX	52	9	43	9	43
Na-ra	IDEXX	9	3	6	3	6

animal hospital						
Hyun-dai animal hospital	Microscopic examination	15	0	15	0	15
Neul-pu-leun animal hospital	Microscopic examination	50	0	50	0	50
Heo animal hospital	Microscopic examination	50	0	50	0	50
Total		176	12	164	12	164

	IDEXX or Microscopic examination			
		Positive	Negative	Total
Anigen Rapid Giardia Ag Test kit	Positive	12	0	12
	Negative	0	164	164
	Total	12	164	176

# (1) Sensitivity

$$\frac{12}{12}$$
 X 100 = 100%

# (2) Specificity

$$\frac{164}{164} \times 100 = 100\%$$