

3-Level Lithography Card Target Value Sheet-20260428

	Unit	Low		Medium		High	
		Standard Value	Upper / Lower Limit (±)	Standard Value	Upper / Lower Limit (±)	Standard Value	Upper / Lower Limit (±)
WBC	10 ⁹ /L	5.98	0.60	9.39	0.94	15.42	1.54
NST	10 ⁹ /L	0.30	0.04	0.47	0.07	0.77	0.12
NSG	10 ⁹ /L	3.47	0.52	5.45	0.82	8.94	1.34
NSH	10 ⁹ /L	0.09	0.02	0.14	0.04	0.23	0.06
LYM	10 ⁹ /L	1.61	0.24	2.54	0.38	4.16	0.62
MON	10 ⁹ /L	0.30	0.04	0.47	0.07	0.77	0.12
EOS	10 ⁹ /L	0.18	0.04	0.28	0.07	0.46	0.12
BAS	10 ⁹ /L	0.03	0.03	0.05	0.05	0.08	0.08
RET	10 ⁹ /L	5.98	0.60	9.39	0.94	15.42	1.54
PLT	10 ⁹ /L	119.60	11.96	187.80	18.78	308.40	30.84
RBC	10 ¹² /L	2.99	0.30	4.70	0.47	7.71	0.77
MCV	fL	75.60	7.56	82.40	8.24	88.60	8.86
HGB	g/L	92.00	9.20	148.00	14.80	191.00	19.10