

Caniv-4

General Information

Purpose : Detection of *Dirofilaria immitis* antigen
Ehrlichia canis antibody
Borrelia burgdorferi antibody
Anaplasma phagocytophilum/*Anaplasma platys* antibody

Principle : Immunochromatographic assay

Specimen : Whole blood, plasma or serum

Materials provided :


- Anigen Rapid Caniv-4 test device
- Assay diluents dropping buffer
- Anticoagulant tube
- Disposable capillary tube

Special Features :

- Simultaneous detection of four kinds of Vector Borne diseases
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RB21-20	Caniv-4	Quattro Device	10 Tests/Kit