

# FHW Ag/Ab

## General Information

**Purpose :** Detection of Feilne heartworm antigen and antibody

**Principle :** Immunochromatographic assay

**Specimen :** Whole blood, plasma or serum

**Materials provided :**


- Anigen Rapid FHW Ag/Ab test device
- Assay diluents dropping bottle
- Anticoagulant tube
- Disposable capillary tube

**Special Features :**

- Simultaneous detection of FHW Ag and Ab with one specimen preparation
- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity



## Test Procedure



## Ordering Information

Cat. No.	Description	Type	Pack size
RC22-05	FHW Ag/Ab	Dual Device	10 Tests/Kit