

Anigen Rapid *New product*

# Total IgE (Allergy indicator test)

## General Information

Allergens drive naïve T cells to differentiate into Th2 cells which can help B cells to undergo heavy chain isotype switching and produce IgE. In other words, IgE signifies that the host might develop hypersensitive immune reaction on subsequent exposure to the allergen. With the Anigen Rapid Total IgE test kit, you can step on to Intradermal Skin Test (IDST)

**Purpose :** Detection of Total IgE

**Principle :** Immunochromatographic assay

**Specimen :** canine serum

**Materials provided :**



- Anigen Rapid Total Ig E test device
- Assay diluents bottle
- Disposable capillary tube
- Disposable dropper

**Special Features :**

- Detection of Immediate Hypersensitivity
- One-step testing procedure : Quick & Accurate
- High Sensitivity, High Specificity



## Test Procedure



### Interpretation

C	T	Negative
C	T	Positive
C	T	Invalid
C	T	Invalid

## Ordering Information

Cat. No.	Description	Type	Pack size
RB2110DD	Total IgE	Device	10 Tests/Kit