

Toxoplasma Ab

General Information

Purpose : Detection of *Toxoplasma gondii* antibody

Principle : Immunochromatographic assay

Specimen : Whole blood, plasma or serum

Materials provided :


- Anigen Rapid Toxo test device
- Assay diluents dropping bottle
- Anticoagulant tube
- Disposable capillary tube

Special Features :

- One-step testing procedure : Quick & Accurate
- High sensitivity and specificity
 - Sensitivity : 99.9% vs. ELISA
 - Specificity : 98.7% vs. ELISA



Test Procedure



Ordering Information

Cat. No.	Description	Type	Pack size
RB22-06	Toxoplasma Ab	Device	10 Tests/Kit