

UC-32B VET

Cute POC Urine Analyzer

Technical Specifications

Parameters

13+1, WBC, KET, URO, BIL, PRO, GLU, SG, BLD, pH, VC, CR, CA, MA + PCR

Principle

Dual Wavelength Reflectance Photometry (470 nm, 550 nm, 620 nm, 720 nm)

Operation System

LCD Screen with Touch Key

Dimension

130*70*29 mm, 200 g

Data Output

USB, Bluetooth

Throughput

Up to 80 T/H

Power Supply

Rechargeable Battery (300 Tests Available per Full Charge)

Reagent Strips*

	WBC	KET	BLD	pH	PRO	GLU	NIT	BIL	SG	URO	VC	MA	CR	Ca
UC VET13 Plus	•	•	•	•	•	•		•	•	•	•	•	•	•
UC VET12	•	•	•	•	•	•	•	•	•	•		•	•	
UC VET11	•	•	•	•	•	•	•	•	•	•	•			
UC VET10	•	•	•	•	•	•	•	•	•	•				


*Strips especially for vet / Customizable strips

URIT Medical Electronic Co., Ltd.

Head Office
Add.: No. D-07 Information Industry District, High-tech Zone,
Guilin, Guangxi 541004, P. R. China
Tel: +86-773-2288586 Fax: +86-773-2288560
E-mail: export@uritest.com service@uritest.com


Shenzhen Office
Add.: Room EJ, 24th Floor, New Baohui Building, Nanhai Road,
Nanshan District, Shenzhen, P.R. China
Tel: +86-755-26413721
http://www.urit.com

URIT Vet Since 2004
All for Vet, All for Better



UC-32B VET
Cute POC Urine Analyzer


Useful, PCR Ratio



Better Diagnosis for Kidney Disease

Efficient, Up to 13+1 Parameters

- | | | |
|-------|-------|-------|
| ■ WBC | ■ VC | ■ GLU |
| ■ KET | ■ PRO | ■ BIL |
| ■ BLD | ■ SG | ■ URO |
| ■ Ca | ■ pH | |
| ■ MA | ■ CR | |
| ☆ PCR | | |




Easy to Use

Big LCD Screen with Touch Key


Bluetooth Connection

Multiple Language Available

Exclusive Software for Veterinary Use



Auto-recognition of Reagent Strips



Put and Test Directly