

WSAVA Vaccination Guideline 2015



WSAVA
Global Veterinary Community





Fig.1 Flow chart for serological testing of puppies



Reference: WSAVA Vaccination Guidelines 2015, p.9 ~ p.10

BioNote's CPV Ab, CDV Ab Rapid Test Kit



According to WSAVA Vaccination Guidelines, dog's protection against several infectious diseases can be diagnosed by rapid antibody test kits. With the antibody test results, vets will be able to set up vaccination schedule and make a prognosis.

BioNote's Rapid CPV Ab, CDV Ab test kits helps you to measure serum antibody titers in a very simple and convenient way.

Meaning of measuring serum anti-CPV Ab, anti-CDV Ab titers

1. Monitor Immunity to canine infectious diseases
2. Detect the presence of protective Ab derived from vaccination (the sample should be taken at least 4 weeks after the last vaccination)
3. Vaccination Scheduling
4. Prevent the spreading of the infectious diseases

