



Veterinary Office Laboratory Testing Introduces the




# Rodenticide Stx™

## The Only 5 Minute Test For Anti-Coagulant Rodenticides

The VOLT Rodenticide Test is a rapid, sensitive test for the semi-quantitative detection of ingested compounds such as Warfarin, Coumadin and all of their derivatives in the anticoagulant class of Rodenticides. This is not a coagulation test, but a self contained test which uses only 75uL of serum or plasma and requires no additional investment

### Specifications:


- 5 minute Go/No Go result
- 75 uL sample without the need for a separate instrument
- Results in nanograms
- LOD is 75 billionths of a gram
- Identification for presence of an Anti-coagulant or derivative
- Screens for Super Rodenticides
- 14 - 16 month expiration



(Shown with optional 75uL reusable pipette)



(Positive & Negative controls available too)



**KACEY VOLT Rodenticide-STX™** are packaged in either 5 or 10 test packs

for your convenience. Also available in our new **Poison Test Station** system with our Ethylene Glycol, Organo-Phosphate tests and related antidotes.

Please visit our web site [www.KaceyDiagnostics.com](http://www.KaceyDiagnostics.com) for full package insert and white papers.

Contact your local KACEY® dealer for more information about the VOLT Rodenticide Stx test.

Authorized Dealer:



**KACEY™**  
DIAGNOSTICS

1854 A Hendersonville Rd  
Asheville, NC 28803

Ph: 1(828)685-3569

Fax: 1(828)685-3571

[www.kaceydiagnostics.com](http://www.kaceydiagnostics.com)

*“Kacey Products Providing Better Animal Care”*